



Fannin Ltd.  
Fannin House, South County Business Park,  
Leopardstown, Dublin 18.  
Tel: + 353 (01) 2907000  
Fax: + 353 (01) 2907118  
[www.fannin.ie](http://www.fannin.ie)



Fannin Ltd.  
I.D.A. Ballybrit Upper Industrial Estate,  
Ballybrit, Galway.  
Tel: + 353 (01) 2907204/5  
Fax: + 353 (091) 771187  
[www.lipdiagnostic.com](http://www.lipdiagnostic.com)

28/06/2017  
LIP Product Information  
W11103 Brilliance UTI Clarity agar

Dear Customer

We would like to advise you that the testing of the above product has been reviewed and a change introduced under our change control process.

QC testing criteria for growth promotion will now comply with manufacturer's recommendations and ISO 12322:1999 for growth of *Proteus mirabilis*.

Inoculation with *Proteus mirabilis* will be NCTC 10975 strain going forward; rather than the previous ATCC 29906. Growth and expected colour reaction of *Proteus mirabilis* is a brown colony with brown halo and no swarming evident.

There has been no change to the product.

Should you require any further information on this product please do not hesitate to contact us.

Marian Quigney  
Quality Control Manager  
00353 (1) 2907210  
[marian.quigney@fannin.eu](mailto:marian.quigney@fannin.eu)

Andre Buckley  
Production Manager,  
00353 (1) 2907206  
[andre.buckley@fannin.eu](mailto:andre.buckley@fannin.eu)