



Fannin Ltd.
 Fannin House, South County Business Park,
 Leopardstown, Dublin 18.
 Tel: + 353 (01) 2907000
 Fax: + 353 (01) 2907118
 www.fannin.ie

Fannin Ltd.
 I.D.A. Industrial Estate,
 Ballybrit Upper, Galway.
 Tel: + 353 (01) 2907204
 www.lipdiagnostic.com

Customer Notification

Update to Antimicrobial Susceptibility Testing

Date of Notification:

20/03/2024

Fannin Ltd., Galway (LIP) Update Announcement

Dear Customer,

Fannin Ltd., Galway (LIP) wish to notify customers of the update to the Antimicrobial Susceptibility Testing panels and certificate of analysis details as outlined below to align with current EUCAST and CLSI guidelines:

Product	W Number	Previous Report Reference	Updated Report Reference
MUELLER HINTON II AGAR 4mm	W11153	Staphylococcus aureus ATCC® 25923	Staphylococcus aureus ATCC® 29213
MUELLER HINTON AGAR (4mm)	W11154	Staphylococcus aureus ATCC® 25923	Staphylococcus aureus ATCC® 29213
MUELLER HINTON + 2% SALT (4mm)	W11155	Staphylococcus aureus ATCC® 25923	Staphylococcus aureus ATCC® 29213
MUELLER HINTON + 5% SHEEP BLOOD (4mm)	W11159	Streptococcus pneumoniae ATCC® 49619 Oxacillin 1 µg specification 8mm-14mm Batch performance tested as per Fannin Ltd standard operational procedures and in accordance with EUCAST methodologies.	Streptococcus pneumoniae ATCC® 49619 Oxacillin 1 µg specification 0mm-12mm Batch performance tested as per Fannin Ltd standard operational procedures and in accordance with CLSI methodologies.
MH AGAR + 5% BLOOD & NAD (4mm) (EUCAST)	W11171	Haemophilus influenzae NCTC® 8468	Haemophilus influenzae ATCC®49766

The above changes will come into effect for all products tested from the 20/03/2024.

This notice needs to be passed on to all those who need to be made aware within your organisation.

Should you have any queries regarding this change please do not hesitate to contact us.

Shane Wilmoth

Shane Wilmoth
Interim QC Manager,
Fannin Ltd, Galway.
Email: shane.wilmoth@fannin.eu

Laura Gallagher

Laura Gallagher
Site Quality Manager,
Fannin Ltd, Galway.
Email: lipquality@fannin.eu